





CD Burning Instructions – Deep Burner (freeware)

(Prerequisite: download the free version of ‘Deep Burner’ CD creation software from www.deepburner.com do not change any configurations when installing simply click continue through all configuration screens)

It should firstly be noted that to burn CD’s from the eRAD system that a CD burning license is required along with the initiating user having forwarding rights and lastly the PC being used must also have a CD drive capable of burning CD’s.

1. Find required study on a worklist.
2. Select the forward function from the tools to the left of the worklist. 
 - multiple studies can be selected and the batch forwarding option should be selected.
3. Select the check box next to ‘Create DICOM CD’.
4. Push the ‘Forward’ button.
5. You will be presented with a box title of ‘Media Options’
 - The first option “Add Viwer” will add a copy of the eRAD diagnostic DICOM viewer to the CD being created
 - The second option “Omit Structured Report objects from DICOMDIR” should only be selected if you do not want any reports associated with the study to be included on the CD
 - The last option is “Uncompress all images” this will copy images in an uncompressed format to allow full image resolution and quality.
6. Press the continue button.
7. An ISO file will now be created. Once this is complete select the “Get ISO file” option, you can then save the file to a given location if using internet explorer or to your desk top if using Mozilla/Firefox. **Remember the name and location of this file. (saving to desk top or a folder on the desk top is a good option)**
8. Select the “Back to list” button.
9. Open the Deep Burner software
10. Select ‘Burn ISO image’ from the ‘Select project type’ prompt
11. Within in the Box labeled “Image location” is the following button: 
select this button and then the location of the ISO file you just saved.
12. Lastly select the “Burn CD” button.